

# INCUBATOR SHAKER



## MODEL SPECIFICATION :

MODEL	SIZE (mm) W X D X H	TEMPERATURE
PIS-1	660 X 765 X 650 mm	3°C Above ambient to 60°C
PIS-2	660 X 765 X 650 mm	5°C to 60°C

## TECHNICAL SPECIFICATION :

Internal body material	Stainless Steel - 304 / 316 Grade
External body material	Powder coated CRCA. Steel (standard models ) / Stainless Steel - 304 grade
Temperature control	Microprocessor with PT-100 sensor
Display	4" LCD screen, large size display for ease of reading
Alarm	Audio visual alarm
Illumination	8 watts fluorescent tube
Insulation (CFC Free polyurethane foam)	70 mm minimum for body & 80mm for door
Noise level	Less than 65 (A)