



Catalogue

PGB220, PGB321
220g, 330g/0.001g


Precision Balance



2024-25

PGB220, 321

CATALOGUE 2024-25

**220g,330g/1mg****Precision Balance****PGB220, PGB321****SALIENT FEATURES**

- ✦ Low Battery Indication
- ✦ Auto Calibration with External Weight
- ✦ Multifunction Weighing Units
- ✦ Auto Zero tracking
- ✦ Overload Protection Design
- ✦ S.S. Weighing Pan, Level indicator, Adjustable Fit
- ✦ High Resolution, Quick Weighing, Accurate Result
- ✦ In built rechargeable battery for continuous use.



E - Catalogue

HSN Code - 90160010

SPECIFICATIONS	MODEL	PGB220	PGB321
	Capacity	220g	330g
	Readability	0.001g (1 mg)	
	Repeatability	0.001g	
	Linearity (±)	0.003g	
	Make	WENSAR	
	Pan Size	90 mm Ø	
	Response Time	3 - 4 Sec	
	Calibration	Automatic External	
	Display	LCD with back light	
	Tare Range	Full	
	Operating Temperature	10°C to 40°C	
	Power supply	Adaptor, Using AC 220V / 50Hz	
	Dimension (LxWxH)	360 x 260 x 355 mm	
	Weight	3.5 kg	

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Adaptor, PAN, Test Certificate.



WENSAR®
AN ISO 9001:2015 COMPANY

Analytical Balance & Standard Weights

WENSAR Products



Semi Micro Balance-0.01mg



Semi Micro Balance-0.01mg



Analytical Balance - 0.1mg



Analytical Balance - 0.1mg




Analytical Balance - 0.1mg



Moisture Analyzer -110g/1mg




Moisture Analyzer -160g/1mg
Internal Calibration



Moisture Analyzer -210g/1mg




Moisture Analyzer -110g/0.1mg



High Precision Balance - 1mg



High Precision Balance -10mg



Precision Balance - 200g/1mg




200g, 300g / 1mg



200g, 500g / 1mg




220g, 320g / 1mg



Density Balance - 220g / 1mg




600g, 1000g / 0.01g



600g / 0.01g



600g / 0.01g



1200g / 0.01g




3000g / 0.01g




3010g / 6000g - 0.1g



3010g / 6000g - 0.1g



5000g / 0.01g



30kg / 0.1g



3kg/0.1g 15kg/0.5g 30kg/1g



20kg / 0.1g



100kg / 1g



60kg/5g,100kg/10g,300kg/20g



Density Kit



Thermal Printer



Extra Display



mg - Wire / Sheet Weights



Standard Weights-Class E1, E2



Standard Weights - F1 Class