


HOT PLATE



MODEL SPECIFICATION : CIRCULAR

MODEL	SIZE (inch)	RATING(Watts)
PHPC-1	8"	1000
PHPC-2	10"	1200
PHPC-3	12"	1500

MODEL SPECIFICATION : RECTANGULAR

MODEL	SIZE (inch)	RATING(Watts)
PHPR-1	10"X12"	1000
PHPR-2	10"X16"	1200
PHPR-3	12"X18"	1500
PHPR-4	18"X24"	2500

TECHNICAL SPECIFICATION :

Plate design	Rectangular & Circular
Cabinet MOC	Powder coated MS
Heating plate surface	Stainless Steel 304 / Mild Steel
Temperature Range	Maximum 300°C, Resolution 1°C
Temperature accuracy	± 1°C
Input voltage	220 V AC ± 10%, 50-60 Hz
Controller	Thermostatic Controller / Digital PID Controller