

equitron



# CONDENSED catalogue



## equitron<sup>®</sup> Water Still - Fully SS



- Industrial grade Incoloy sheathed immersion heater
- Dry run safety alarm & cut off
- Can be wall mounted or table top
- Supplied complete with raw water flow rate meter and tubings
- With spare gasket



Model	Inner Dimensions	Capacity
#7655M	Water Still with input water flow meter & baffle	4 L / hr

## lapiz<sup>®</sup> Colony Counter



- Freedom from wired pens, use any commercially available hard tipped marker pen
- Circular tube for excellent all round illumination
- Powerful magnifying lens reduces operator fatigue
- Fully stainless steel GMP model
- With last count error decrement switch



Model	Product	Count upto
#0671M	Digital - desk model - SS GMP model	9,999

## equitron<sup>®</sup> Slide Cabinet – Individual



- Expandable & Flexible storage system for slides
- Each stackable cabinet stores 2520 slides (individual) in just 1.25<sup>2</sup> ft
- Can stack upto 3 cabinets - one above the other
- Drawer does not fall down even when fully pulled out



Model	Product	Capacity
#6005	Slide Cabinet - Individual Storage - 36 drawers	2520 slides
#6050	Base for above	Stack upto 3 cabinets

## equitron<sup>®</sup> Slide Staining Racks & Troughs



- Stainless steel die punched rack
- Stainless steel joint less staining trough with rounded edges for easy cleaning



Model	Product	Capacity
#6142	Slide Staining Rack - SS	10 slides
#6145	Slide Staining Rack - SS	25 slides
#6152	Slide Staining Trough - SS	10 slides
#6155	Slide Staining Trough - SS	25 slides

## equitron<sup>®</sup> Slide Trays / Stand



- Viewing tray: aluminium embossed and powder coated
- Drying tray: stainless steel glossy finish - drains off water easily
- Staining Stand (ghodi): full stainless steel rigid construction



Model	Product	Capacity / Description
#6105	Viewing tray - AL	24 slides
#6115	Drying tray - SS	4 grooves of 300 mm
#6121	Staining Stand (Ghodi) - SS	280 mm long