





MOISTURE ANALYZER

PGB 1MB


Salient Features

- Halogen Light Heating
- Curve (Graph)
- Temperature and time can be set
- The percentage of Moisture content
- The percentage of dry residual
- Touch Screen
- Store historical setting
- Loadcell Sensor



MOISTURE DETERMINATION (1)

Machine Type : PGB 1MB
Heating Mode : STANDARD
Stop Mode : AUTO STOP
Heating Temp. : 130°C
Elasapsed Time : 03:25
MW : 2.056g
DW : 1.933g
Moisture : 5.98%M
Date : 2017-10-17
Time : 08:16:14



STANDARD ACCESSORIES

The sample pan holder,
Tray handle,
Power cord,
Aluminum sample pan
Instruction Manual

Moisture Analyzer

SPECIFICATIONS	MODEL	PGB 1MB
	Capacity	110g
	Make	WENSAR
	Readability	0.001 g (1mg)
	Repeatability (±)	0.001 g
	Linearity (±)	0.002 g
	Display	Touch Screen
	Pan Size	90 mm Ø
	Sensor	loadcell
	Calibration	Automatic External
	Operating Temperature	5°C to 35°C
	Heating Source	Halogen lamp
	Temp. Sensor	PT - 100
	Temperature Setting	40°C~199°C
	Temp. Step	1°C
	Time Setting	1~99 minutes
	Dry Readability	0.01%
	Dry Range	100.00% ~ 0.00%
	Curve	Yes
	Allumium Sample PAN	Available
	No. of Storage	20
	Power Supply	AC 220V / 50-60 Hz
	Dimension (LxWxH) mm	490 x 360 x 350 mm
	Weight	7.5 kg

★ Specifications are subject to change without notice.